


Sterling™

PTA Balloon Dilatation Catheter



Demand

from your 0.018" / 0.46 mm Balloon Dilatation Catheter

More

The #1 0.018"/ 0.46 mm**
PTA Balloon Dilatation Catheter

with the most comprehensive
size offering available

2 - 10 mm diameters
10 - 220 mm lengths*

The most comprehensive
size offering available

10 seconds*

Best-in-class
deflation time

0.020"/
0.51 mm*

Lowest in-class
lesion entry profile

The Sterling™ family† is designed to offer best-in-class lesion entry profile and deflation times with the most clinically comprehensive size offering and indications available.



Demand Exceptional Crossability

- Exceptionally low 0.020"/ 0.51 mm lesion entry profile on a 0.018"/ 0.46 mm wire
- 7x220 through 5F / 0.65 mm



Demand the Broadest Size Matrix

- The most comprehensive size offering available in both Over-The-Wire and Monorail™ rapid-exchange
- Indicated for wider acute use including carotid‡ and post stent dilatation

* Average measurements taken by Boston Scientific (n=3, 5 x 150 mm balloons). Data on file.
Bench testing may not be representative of clinical performance.

** US and EU market share data per internal and external analysis.

† The Sterling family includes Sterling and Sterling SL.

‡ Monorail only

Sterling™ and Sterling™ SL Over-The-Wire Balloon Dilatation Catheter

Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
2.0	20	90	4	H74939032200290	
2.0	20	150	4	H74939032200210	
2.0	30	90	4	H74939032200390	
2.0	30	150	4	H74939032200310	
2.0	40	90	4	H74939032200490	
2.0	40	150	4	H74939032200410	
2.0	60	90	4	H74939032200690	
2.0	60	150	4	H74939032200610	
2.0	80	90	4	H74939148208090**	
2.0	80	150	4	H74939148208010**	
2.0	100	90	4	H74939148201090**	
2.0	100	150	4	H74939148201010**	
2.0	120	90	4	H74939148201290**	
2.0	120	150	4	H74939148201210**	
2.0	150	90	4	H74939148201590**	
2.0	150	150	4	H74939148201510**	
2.0	220	90	4	H74939032202290	
2.0	220	150	4	H74939032202210	
2.5	20	90	4	H74939032250290	
2.5	20	150	4	H74939032250210	
2.5	30	90	4	H74939032250390	
2.5	30	150	4	H74939032250310	
2.5	40	90	4	H74939032250490	
2.5	40	150	4	H74939032250410	
2.5	60	90	4	H74939032250690	
2.5	60	150	4	H74939032250610	
2.5	80	90	4	H74939148258090	
2.5	80	150	4	H74939148258010**	
2.5	100	90	4	H74939148251090**	
2.5	100	150	4	H74939148251010**	
2.5	120	90	4	H74939148251290**	
2.5	120	150	4	H74939148251210**	
2.5	150	90	4	H74939148251590**	
2.5	150	150	4	H74939148251510**	
2.5	220	90	4	H74939032252290	
2.5	220	150	4	H74939032252210	
3.0	20	90	4	H74939032300290	
3.0	20	150	4	H74939032300210	
3.0	30	90	4	H74939032300390	
3.0	30	150	4	H74939032300310	
3.0	40	90	4	H74939032300490	
3.0	40	150	4	H74939032300410	
3.0	60	90	4	H74939032300690	
3.0	60	150	4	H74939032300610	
3.0	80	90	4	H74939148308090**	
3.0	80	150	4	H74939148308010**	
3.0	100	90	4	H74939148301090**	
3.0	100	150	4	H74939148301010**	

Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
3.0	120	90	4	H74939148301290**	
3.0	120	150	4	H74939148301210**	
3.0	150	90	4	H74939148301590**	
3.0	150	150	4	H74939148301510**	
3.0	220	90	4	H74939032302290	
3.0	220	150	4	H74939032302210	
3.5	20	90	4	H74939032350290	
3.5	20	150	4	H74939032350210	
3.5	30	90	4	H74939032350390	
3.5	30	150	4	H74939032350310	
3.5	40	90	4	H74939032350490	
3.5	40	150	4	H74939032350410	
3.5	60	90	4	H74939032350690	
3.5	60	150	4	H74939032350610	
3.5	80	90	4	H74939148358090**	
3.5	80	150	4	H74939148358010**	
3.5	100	90	4	H74939148351090**	
3.5	100	150	4	H74939148351010**	
3.5	120	90	4	H74939148351290**	
3.5	120	150	4	H74939148351210**	
3.5	150	90	4	H74939148351590**	
3.5	150	150	4	H74939148351510**	
3.5	220	90	4	H74939032352290	
3.5	220	150	4	H74939032352210	
4.0	20	40	4	H74939032402040	
4.0	20	80	4	H74939032402080	
4.0	20	135	4	H74939032402010	
4.0	30	80	4	H74939032403080	
4.0	30	135	4	H74939032403010	
4.0	40	40	4	H74939032404040	
4.0	40	80	4	H74939032404080	
4.0	40	135	4	H74939032404010	
4.0	60	80	4	H74939032406080	
4.0	60	135	4	H74939032406010	
4.0	80	80	4	H74939032408080	
4.0	80	135	4	H74939032408010	
4.0	80	150	4	H74939148408010**	
4.0	100	80	4	H74939032401080	
4.0	100	135	4	H74939032401010	
4.0	100	150	4	H74939148401010**	
4.0	120	90	4	H74939148401290**	
4.0	120	150	4	H74939148401210**	
4.0	150	90	4	H74939148401590**	
4.0	150	150	4	H74939148401510**	
4.0	220	90	4	H74939032402290	
4.0	220	150	4	H74939032402210	
5.0	20	40	4	H74939032502040	
5.0	20	80	4	H74939032502080	

Sterling™ and Sterling™ SL Over-the-Wire Balloon Dilatation Catheter (continued)

Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
5.0	20	135	4	H74939032502010	
5.0	30	80	4	H74939032503080	
5.0	30	135	4	H74939032503010	
5.0	40	40	4	H74939032504040	
5.0	40	80	4	H74939032504080	
5.0	40	135	4	H74939032504010	
5.0	60	80	4	H74939032506080	
5.0	60	135	4	H74939032506010	
5.0	80	80	4	H74939032508080	
5.0	80	135	4	H74939032508010	
5.0	100	80	4	H74939032501080	
5.0	100	135	4	H74939032501010	
5.0	120	90	5	H74939032501290	
5.0	120	150	5	H74939032501210	
5.0	150	90	5	H74939032501590	
5.0	150	150	5	H74939032501510	
5.0	200	90	5	H74939032502090	
5.0	200	150	5	H74939032502110	
5.0	220	90	5	H74939032502290	
5.0	220	150	5	H74939032502210	
6.0	20	40	4	H74939032602040	
6.0	20	80	4	H74939032602080	
6.0	20	135	4	H74939032602010	
6.0	30	80	4	H74939032603080	
6.0	30	135	4	H74939032603010	
6.0	40	40	4	H74939032604040	
6.0	40	80	4	H74939032604080	
6.0	40	135	4	H74939032604010	
6.0	60	80	4	H74939032606080	
6.0	60	135	4	H74939032606010	
6.0	80	80	5	H74939032608080	
6.0	80	135	5	H74939032608010	
6.0	100	80	5	H74939032601080	
6.0	100	135	5	H74939032601010	
6.0	120	90	5	H74939032601290	
6.0	120	150	5	H74939032601210	
6.0	150	90	5	H74939032601590	
6.0	150	150	5	H74939032601510	
6.0	200	90	5	H74939032602090	
6.0	200	150	5	H74939032602110	
6.0	220	90	5	H74939032602290	
6.0	220	150	5	H74939032602210	
7.0	20	40	4	H74939032702040	
7.0	20	80	4	H74939032702080	
7.0	20	135	4	H74939032702010	
7.0	30	80	4	H74939032703080	
7.0	30	135	4	H74939032703010	
7.0	40	40	4	H74939032704040	

Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
7.0	40	80	4	H74939032704080	
7.0	40	135	4	H74939032704010	
7.0	60	80	5	H74939032706080	
7.0	60	135	5	H74939032706010	
7.0	80	80	6	H74939032708080	
7.0	80	135	6	H74939032708010	
7.0	100	80	6	H74939032701080	
7.0	100	135	6	H74939032701010	
7.0	120	90	5	H74939032701290	
7.0	120	150	5	H74939032701210	
7.0	150	90	5	H74939032701590	
7.0	150	150	5	H74939032701510	
7.0	200	90	5	H74939032702090	
7.0	200	150	5	H74939032702110	
7.0	220	90	5	H74939032702290	
7.0	220	150	5	H74939032702210	
8.0	20	40	4	H74939032802040	
8.0	20	80	4	H74939032802080	
8.0	20	135	4	H74939032802010	
8.0	30	80	4	H74939032803080	
8.0	30	135	4	H74939032803010	
8.0	40	40	4	H74939032804040	
8.0	40	80	4	H74939032804080	
8.0	40	135	4	H74939032804010	
8.0	60	80	5	H74939032806080	
8.0	60	135	5	H74939032806010	
8.0	80	80	6	H74939032808080	
8.0	80	135	5	H74939032808010	
9.0	20	80	5	H74939032902080	
9.0	20	135	5	H74939032902010	
9.0	30	80	5	H74939032903080	
9.0	30	135	5	H74939032903010	
9.0	40	80	5	H74939032904080	
9.0	40	135	5	H74939032904010	
9.0	60	80	5	H74939032906080	
9.0	60	135	5	H74939032906010	
9.0	80	80	6	H74939032908080	
9.0	80	135	6	H74939032908010	
10.0	20	80	5	H74939032102080	
10.0	20	135	5	H74939032102010	
10.0	30	80	5	H74939032103080	
10.0	30	135	5	H74939032103010	
10.0	40	80	5	H74939032104080	
10.0	40	135	5	H74939032104010	
10.0	60	80	6	H74939032106080	
10.0	60	135	6	H74939032106010	
10.0	80	80	6	H74939032108080	
10.0	80	135	6	H74939032108010	

* 4F = 1.33mm, 5F = 1.7mm, 6F = 2.0mm

**Sterling™ SL

Sterling™ and Sterling™ SL Monorail Balloon Dilatation Catheter

Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
2.0	20	150	4	H74939031202010	
2.0	40	150	4	H74939031204010	
2.0	60	150	4	H74939031206010	
2.0	80	90	4	H74939147208090**	
2.0	80	150	4	H74939147208010**	
2.0	100	90	4	H74939147201090**	
2.0	100	150	4	H74939147201010**	
2.0	120	90	4	H74939147201290**	
2.0	120	150	4	H74939147201210**	
2.0	150	90	4	H74939147201590**	
2.0	150	150	4	H74939147201510**	
2.0	220	150	4	H74939031202210	
2.5	20	150	4	H74939031252010	
2.5	40	150	4	H74939031254010	
2.5	60	150	4	H74939031256010	
2.5	80	90	4	H74939147258090**	
2.5	80	150	4	H74939147258010**	
2.5	100	90	4	H74939147251090**	
2.5	100	150	4	H74939147251010**	
2.5	120	90	4	H74939147251290**	
2.5	120	150	4	H74939147251210**	
2.5	150	90	4	H74939147251590**	
2.5	150	150	4	H74939147251510**	
2.5	220	150	4	H74939031252210	
3.0	20	135	4	H74939031302010	
3.0	30	135	4	H74939031303010	
3.0	40	135	4	H74939031304010	
3.0	60	135	4	H74939031306010	
3.0	80	90	4	H74939147308090	
3.0	80	150	4	H74939147308010**	
3.0	100	90	4	H74939147301090**	
3.0	100	150	4	H74939147301010**	
3.0	120	90	4	H74939147301290**	
3.0	120	150	4	H74939147301210**	
3.0	150	90	4	H74939147301590**	
3.0	150	150	4	H74939147301510**	
3.0	220	150	4	H74939031302210	
3.5	20	135	4	H74939031352010	
3.5	30	135	4	H74939031353010	
3.5	40	135	4	H74939031354010	
3.5	60	150	4	H74939031356010	
3.5	80	90	4	H74939147358090**	
3.5	80	150	4	H74939147358010**	
3.5	100	90	4	H74939147351090**	
3.5	100	150	4	H74939147351010**	
3.5	120	90	4	H74939147351290**	
3.5	120	150	4	H74939147351210**	
3.5	150	90	4	H74939147351590**	
3.5	150	150	4	H74939147351510**	
3.5	220	150	4	H74939031352210	
4.0	10	80	4	H74939031401080	
4.0	15	80	4	H74939031401580	
4.0	20	80	4	H74939031402080	
4.0	20	135	4	H74939031402010	
4.0	30	80	4	H74939031403080	
4.0	30	135	4	H74939031403010	
4.0	40	80	4	H74939031404080	
4.0	40	135	4	H74939031404010	
4.0	60	135	4	H74939031406010	
4.0	80	90	4	H74939147408090**	
4.0	80	150	4	H74939147408010**	
4.0	100	90	4	H74939147401090**	
4.0	100	150	4	H74939147401010**	


Diameter (mm)	Length (mm)	Shaft (cm)	Sheath Compatibility (F)*	UPN Number	Qty
4.0	120	90	4	H74939147401290**	
4.0	120	150	4	H74939147401210**	
4.0	150	90	4	H74939147401590**	
4.0	150	150	4	H74939147401510**	
4.0	220	150	4	H74939031402210	
4.5	15	80	4	H74939031451580	
4.5	20	80	4	H74939031452080	
4.5	20	135	4	H74939031452010	
4.5	30	80	4	H74939031453080	
4.5	30	135	4	H74939031453010	
4.5	40	80	4	H74939031454080	
4.5	40	135	4	H74939031454010	
5.0	10	80	4	H74939031501080	
5.0	15	80	4	H74939031501580	
5.0	20	80	4	H74939031502080	
5.0	20	135	4	H74939031502010	
5.0	30	80	4	H74939031503080	
5.0	30	135	4	H74939031503010	
5.0	40	80	4	H74939031504080	
5.0	40	135	4	H74939031504010	
5.0	60	135	4	H74939031506010	
5.5	15	80	4	H74939031551580	
5.5	20	80	4	H74939031552080	
5.5	20	135	4	H74939031552010	
5.5	30	80	4	H74939031553080	
5.5	30	135	4	H74939031553010	
5.5	40	80	4	H74939031554080	
5.5	40	135	4	H74939031554010	
6.0	10	80	4	H74939031601080	
6.0	15	80	4	H74939031601580	
6.0	20	80	4	H74939031602080	
6.0	20	135	4	H74939031602010	
6.0	30	80	4	H74939031603080	
6.0	30	135	4	H74939031603010	
6.0	40	80	4	H74939031604080	
6.0	40	135	4	H74939031604010	
6.0	60	135	4	H74939031606010	
6.5	15	80	4	H74939031651580	
6.5	20	80	4	H74939031652080	
6.5	20	135	4	H74939031652010	
6.5	30	80	4	H74939031653080	
6.5	30	135	4	H74939031653010	
6.5	40	80	4	H74939031654080	
6.5	40	135	4	H74939031654010	
7.0	10	80	4	H74939031701080	
7.0	15	80	4	H74939031701580	
7.0	20	80	4	H74939031702080	
7.0	20	135	4	H74939031702010	
7.0	30	80	4	H74939031703080	
7.0	30	135	4	H74939031703010	
7.0	40	80	4	H74939031704080	
7.0	40	135	4	H74939031704010	
7.0	60	135	5	H74939031706010	
8.0	10	80	4	H74939031801080	
8.0	15	80	4	H74939031801580	
8.0	20	80	4	H74939031802080	
8.0	20	135	4	H74939031802010	
8.0	30	80	4	H74939031803080	
8.0	30	135	4	H74939031803010	
8.0	40	80	4	H74939031804080	
8.0	40	135	4	H74939031804010	
8.0	60	135	5	H74939031806010	

* 4F = 1.33 mm, 5F = 1.7 mm, 6F = 2.0 mm

** Sterling™ SL


Sterling™ Monorail™

PTA Balloon Dilatation Catheter



Sterling™ Over-the-Wire

PTA Balloon Dilatation Catheter



		Length (mm)											
Diameter (mm)	Balloon	10	15	20	30	40	60	80	100	120	150	200	220
	2.0 mm												
	2.5 mm												
	3.0 mm												
	3.5 mm												
	4.0 mm												
	4.5 mm												
	5.0 mm												
	5.5 mm												
	6.0 mm												
	6.5 mm												
	7.0 mm												
	8.0 mm												
	9.0 mm												
	10.0 mm												

		Shaft Length (cm)			
Platform		40	80/90	135	150
OTW					
Rapid Exchange					

Boston Scientific

Advancing science for life™

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

PI-273018-AA / PI-195601-AB. Printed in Germany by medicalvision.